

1 inch



- Full Face Delivery
- Direct Nasal Delivery
- Nasal Delivery
- Oral Delivery



Infinity⁴⁸¹
DIRECT NASAL MASK



FlexiFit⁴⁰⁷
PREMIUM NASAL MASK



FlexiFit⁴⁰⁶
PETITE NASAL MASK



CPAPMAN.COM
BEST SERVICE SINCE 1997

Click Here To Order Today

Fisher & Paykel Multi-Mask Sizing Gauge

1. Print document - Note do not choose to scale document to paper size when printing. (Square at top of page should be one inch wide/tall)
2. Cut out paper inside of mask outlines
3. Choose correct size according to fit (406/407) or size lines (431/432/481).

FlexiFit^{431/432}
FULL FACE MASK

185044349REV A

SleepStyle[™]
SLEEP PERFORMANCE FOR
AN ENERGIZED LIFESTYLE

Fisher & Paykel
HEALTHCARE