



**CPAPMAN.COM**  
BEST SERVICE SINCE 1997

***Click Here To Order Today***

**PHILIPS**

Philips Respironics  
1740 Golden Mile Highway  
Monroeville, PA 15146

August 23, 2016

Re: Declaration for PAP and Respiratory Lithium-Ion Battery (PN1126946)

To Whom It May Concern:

We hereby certify that the Lithium-ion Batteries used in the Respironics PAP and Respiratory Lithium ion Battery comply with the following requirement identified in the IATA Lithium Battery Guidance Document which is based on the provisions set out in the 2015-2016 Edition of ICAO Technical Instruction for the Safe Transport of Dangerous Goods by Air and the 56th Edition of the IATA Dangerous Goods Regulations (DGR):

**2.3.4. 7 Lithium Battery-Powered Electronic Devices**

Lithium battery-powered electronic devices are permitted in checked and carry-on baggage with the approval of the operator as follows:

(a) Portable medical electronic devices (Automated External Defibrillators (AED), Nebulizer, Continuous Positive Airway Pressure (CPAP), etc.) containing lithium metal or lithium ion cells or batteries may be carried by passengers for medical use as follows:

3. Batteries must be of a type that meets the requirements of the UN Manual of Tests and Criteria, Part III subsection 38.3.

The batteries have successfully passed the test procedures of Part III, Sub- Section 38.3 of the UN Manual of Tests and Criteria.

Some airlines may consider a maximum wattage for lithium ion batteries. The Respironics PAP and Respiratory lithium ion battery has a Watt Hour rating of approximately 90 Wh.

If the battery pack will be used for airline travel, please consult the airline restrictions, (e.g., number of spare batteries permitted). Please familiarize yourself with airline rules. Before you fly, contact all of your airlines on each of your connecting flights for clarification as to what can be transported via Carry-on and Checked Baggage. Airlines have different policies and it is recommended to request permission to carry your Sleep or Respiratory Care System on board.

In case of incident, accident and/or damage please follow the steps listed in the Material Safety Data Sheet for Transportation, First Aid, General Safety Recommendations and Fire Safety.

For more information, please refer to either the User Instruction Manual for more detailed information regarding the set-up, performance parameters, technical specifications and maintenance of the device and batteries.

Please feel free to contact customer service at 1-800-345-6443 for more information regarding the PAP or Respiratory Lithium ion batteries.

Sincerely,

Philips Respironics