



F&P Evora™ Nasal



F&P Evora Full

		Nose Width		Size	Nose Length	Nose Width
		44 mm or less	More than 44 mm	XS	18 mm or less	Any
Nose Length	17 mm or less	S	W	S-M	More than 18 mm	40 mm or less
	More than 17 mm	M	L	L	More than 18 mm	More than 40 mm



F&P myMask App*

For the full setup experience, download the **F&P myMask™** App or scan the QR code below.



* Not available in all countries.

F&P, Evora and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. 620939 REV A 2021-02 © 2021 Fisher & Paykel Healthcare Limited.

This material is not intended to replace the Evora Nasal or Evora Full sizing guides.