



**CPAPMAN.COM**  
BEST SERVICE SINCE 1997

***Click Here To Order Today***

# *Freedom to move.*

## **Opus<sup>360</sup>**

**NASAL PILLOWS MASK**

The Opus<sup>360</sup> Nasal Pillows Mask is small, quiet and superbly comfortable, offering unprecedented freedom. The slim-line design allows for almost any sleeping position without mask displacement. Its innovative 360-degree pivoting ball-and-socket elbow provides greater versatility. With its new Silicone Pillows and facially Contoured Mask Frame for extra stability, patient comfort is paramount.


Elegant in design and simplicity, let the Opus<sup>360</sup> Nasal Pillows Mask from Fisher & Paykel Healthcare work for you today!

**Contoured Mask Frame  
& "Off-the-Lip" Pillows**



**Fisher & Paykel**  
HEALTHCARE

# Direct Nasal Delivery – Opus<sup>360</sup>




The Opus<sup>360</sup> Nasal Pillows Mask is small, quiet and superbly comfortable, with anatomical pillow-style Silicone Seals and a facially Contoured Mask Frame. A ball-and-socket elbow combines with an innovative anchoring system to give unprecedented freedom of movement and stability. The Opus<sup>360</sup> Nasal Pillows Mask is a complete, Out-of-Box Solution, fully assembled and ready for use with alternative sizes of Silicone Pillows included.

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### Anatomically Shaped Pillows Mask comfort

The Anatomically Shaped Pillows provide easy fitting with an excellent seal. The contours of the nasal pillows ensure greater comfort. The range of sizes provided (S, M, L) ensures comfort and seal are not compromised.




### Contoured Mask Frame Mask comfort and stability

The facially Contoured Mask Frame, with smooth soft pads, conforms naturally to the shape of the face to provide a comfortable fit and reliable support. Its slim-line design allows side-sleeping without mask displacement.




### Tube Anchoring System Mask stability

The innovative pivot system allows the tube to be anchored at the top or side, or left free to track. Elegant in design and simplicity, it firmly holds the tube in the desired position. Alternatively, the tube can be left unanchored so that the ball and socket elbow allows true 360-degree freedom.



### Efficient Exhaust Quiet usage

The advanced air diffuser system is quiet, minimizing the noise and draft of the exhaust flow to prevent both sleeper and partner disturbance.




### "Off-the-lip" Pillows Mask comfort

The Silicone Seal Pillows are soft on the nose and designed to be "off the upper lip" for maximum comfort.



### Premium Frosted Silicone Seal Optimizing seal technology

The enhanced seal technology has a soft, smooth feel and slides easily into position. Its frosted texture optimizes mask seal performance.



### Hidden Strap Release One-step removal

The easy one-step removal of headgear for cleaning enables reassembly without having to make readjustments. Simply slide the Side Strap down and out of the Hidden Strap Release hooks.

## CONTENTS

- Complete Mask with Mask Frame assembled and ready for use
- All three sizes of Silicone Seal included (S, M, L)
- User-friendly Fitting Guide

REORDER CODE: HC482

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

**Manufacturer** Fisher & Paykel Healthcare Limited, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand  
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: [info@fphcare.com](mailto:info@fphcare.com) Web: [www.fphcare.com](http://www.fphcare.com)

**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Northern Ireland** Tel: 0800 132 189 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA** Tel: 1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001