


Note: For proper mask sizing print out this page at 100%.

Do NOT select the "Shrink oversized pages to paper size" option.

The box in the upper right corner should measure 1 inch square.



A large graphic showing four nasal mask templates labeled P, S, M, and L. Each template has a horizontal double-headed arrow across its base with the text "cut out along lines" below it. To the left of the templates are two small diagrams: one showing a person's face with a vertical line for height (L) and a horizontal line for width (W), and another showing a mask with similar L and W dimensions. Below the templates, the text "ComfortGel nasal mask sizing gauge" is centered, with the website "www.philips.com/comfortgel" underneath. The Philips Resironics logo is in the bottom right corner of the gauge area.

The logo for CPAPMAN.COM features a stylized blue figure inside a white outline of a mask. To the right of the logo, the text "CPAPMAN.COM" is written in large blue letters, with "BEST SERVICE SINCE 1997" in smaller orange letters below it. Underneath, the phrase "Click Here To Order Today" is written in red, italicized font.