

## SleepWeaver Élan Sizing Guide

SleepWeaver Élan Soft Cloth Nasal CPAP Masks come in three sizes. The primary difference is the size of the nasal port in the Zzzephyr Seal as shown below. Most adults will use the Regular. If a Regular cushion does not fit, choosing between a Small or Large mask should be fairly easy based on the prominence of the user's nasal bridge or width of nostrils.



**CPAPMAN.COM**  
BEST SERVICE SINCE 1997

[Click Here To Order Today](#)

### Zzzephyr Seal


### Size Identification



Large



Regular



Small

