

# TAP<sup>®</sup> PAP Nasal Pillow Mask

## GET TO KNOW THE TAP<sup>®</sup> PAP NASAL PILLOW MASK

The TAP PAP Nasal Pillow Mask is intended to be used by individuals who have been diagnosed by a Physician as requiring CPAP or Bi-Level ventilator treatment. The TAP PAP Nasal Pillow Mask is intended for single patient adult (>66 lb/30 kg) use in the home, hospital or other clinical setting.



**CPAPMAN.COM**  
BEST SERVICE SINCE 1997


*Click Here To Order Today*

## TECHNICAL SPECIFICATIONS

### BIAS FLOW DATA

Mask Pressure (cm H <sub>2</sub> O)	4	6	8	10	12	14	16	18	20
Vent Flow Rate (L/min)	19	23	27	30	33	35	38	40	42

### BIAS FLOW GRAPH



### RESISTANCE TO FLOW DATA

Drop in pressure (through mask at given flow rate)

- At 50 L/min 0.2 cm H<sub>2</sub>O
- At 100 L/min 0.6 cm H<sub>2</sub>O

Drop in pressure with Medium Seal in place

- At 50 L/min 0.6 cm H<sub>2</sub>O
- At 100 L/min 2.1 cm H<sub>2</sub>O

### MASK DEAD SPACE

- 95 mL (the physical dead space is the empty volume of the mask to the end of the swivel)

### SOUND


- The Sound Power Level of the mask is 25.7 dBA with uncertainty 1.96 dBA.
- The Sound Pressure Level of the mask is 17.8 dBA with uncertainty 1.96 dBA.

## CHINSTRAP




Patients who mouth breathe and mouth leak will benefit from using the Airway Management Chinstrap.

## SPARE PARTS




TAP PAP - Nasal Pillow Mask  
**PAPNP1000**




TAP PAP - Nasal Pillow Mask - No Headgear  
**PAPNP1001**



TAP PAP - NP - Swivel Connector  
**PAPNP1502**




TAP PAP - Nasal Pillow Seal  
Small: **PAPNP1101**  
Medium: **PAPNP1102**  
Large: **PAPNP1103**




TAP PAP - Upper Mouth Piece  
**PAPNP1201**



TAP PAP - Nasal Pillow - Headgear  
**PAPNP1300**



TAP PAP Chin Strap  
**PAP000400**



TAP PAP - Nasal Pillow - Plastic Nut and Bolt  
**PAPNP1500**

To visit our website with your smartphone, open your QR Code Reader App and scan:



3418 Midcourt Road, Suite 114  
Carrollton, TX, 75006, USA

**866.264.7667**