



Click Here To Order Today

At home or on-the-go,
Transcend provides the
sleep apnea therapy you need.

Featuring its small, light and innovative design, Transcend is the gold standard in portable PAP innovation. With a full line of PAP products and accessories unlike any other CPAP on the market, you can customize your Transcend device to fit your active lifestyle. Sleep apnea therapy is no longer limited to the night stand.



SLEEP APNEA THERAPY

Customized therapy
is as easy as ①-②-③

Transcend's patient-centric product line features accessory options unlike any other sleep apnea therapy system on the market to give you therapy that fits your lifestyle.

① Therapy

- Transcend Auto™ with expiratory relief
- Transcend EZEX™ CPAP with expiratory relief
- Transcend® CPAP
- TranSync™ wireless therapy monitoring




② Humidification

- Transcend Heated Humidifier™



③ Power Options

- P8™ Battery*
- P4™ Battery*
- Transcend Portable Solar Battery Charger™
- DC Mobile Power Adaptor



* Not for use with Transcend Heated Humidifier.

**Contact your local Transcend dealer
for more information.**



www.MyTranscend.com © 2013 Somnics International, Inc.
Protected by U.S. and foreign patents issued and pending. Transcend® is a registered trademark of Somnics International, Inc. - PN: 400130 Rev. B

SLEEP APNEA THERAPY



Sleep apnea therapy
for the active lifestyle.



Transcend is changing the way sleep apnea therapy is delivered.™

Identical in size, each Transcend device is designed to fit your active lifestyle.



Transcend AUTO™

Transcend Auto is a fully-featured APAP that continually monitors your breathing to adapt to your changing therapy needs. Transcend EZEX technology provides pressure relief on exhalation for more breathing comfort.



Transcend EZEX™

Transcend EZEX is a fully featured fixed-pressure CPAP with EZEX exhalation pressure relief technology for added comfort.



Transcend®

Transcend is a fully featured fixed-pressure CPAP. This original CPAP is the pioneer in portable sleep apnea therapy. At home or on the road, Transcend is designed to enhance the lifestyles of today's active CPAP users.

Small size, big features for therapy anywhere.

Transcend's fully-featured therapy devices are small, lightweight and portable. They are perfect for packing into a suitcase or taking up minimal space on the night stand.



Items shown proportionate.

	Transcend AUTO™	Transcend EZEX™	Transcend®
Auto-titrating	✓		
Fixed pressure	✓	✓	✓
Operating pressure (4-20 cmH2O)	✓	✓	✓
Transcend EZEX™ pressure relief (settings 1,2,3)	✓	✓	
Ramp (0-45 min), adjustable	✓	✓	✓
Weighs less than one pound (426 g)	✓	✓	✓
Automatic altitude adjustment	✓	✓	✓
AHI & leak detect	✓	✓	✓
Up to 13 months compliance data storage	✓	✓	✓
TranSync™ Wireless Compliance capability	✓	✓	
Compatible with Transcend batteries	✓	✓	✓
Compatible with Transcend Heated Humidifier™	✓	✓	✓

For full specifications, please visit www.myTranscend.com

Contact your medical supplier for more information.