


[Click Here To Order Today](#)

Mirage Vista Mask Fitting Template

How to use the fitting template:

1. Cut out template along the dotted line
2. Follow the directions on the template to determine appropriate cushion size



The diagram shows a large rectangular template with a dashed border. On the right side, there is a pair of scissors icon. The template is divided into two main sections. The top section contains two mask cushion shapes, one labeled 'DEEP' and one labeled 'STANDARD'. Each shape has an upward-pointing arrow and the text 'ALIGN WITH BOTTOM OF NOSE'. The bottom section is divided into two columns. The left column contains the text 'RESMED Mirage Vista™ NASAL MASK Fitting Template' and a note. The right column contains a diagram of a person's face with a mask, and two numbered steps: '1 Check the nostril width for clearance' and '2 Check the nasal bridge width for clearance.' The diagram shows arrows indicating the measurement points for both steps.

RESMED

Mirage Vista™
NASAL MASK


Fitting Template

Note: Mask fitting templates should be used to aide patients in the fitting of CPAP masks. Altering the size or dimensions of the original online PDF in any way may compromise the quality and accuracy of this template. The use of this template is not intended to replace the advice of a medical professional, sleep technician or HME dealer.

Follow these two steps to choose the correct size:

1 Check the nostril width for clearance

2 Check the nasal bridge width for clearance.



Actual Size: 9 1/2 inch (Width) X 4 3/4 inch (Height)

Ruler is for scaling purposes. To ensure this print is accurate, match a ruler with the ruler above.