



CPAPMAN.COM
BEST SERVICE SINCE 1997



ZOPEC MEDICAL™
Innovation for Health

Sleep Safe. Sleep Uninterrupted.

**Explore Travel Batteries
& Uninterruptible Power Supplies**



- + **Highest class of safety** - engineered for hospitals but made for you at home.
- + **Sleep Uninterrupted** - built to be a daily back power supply for uninterrupted therapy.
- + **Certified Globally** - FDA, UL, CE, CSA, IEC, FCC, RoHS.
- + **Guaranteed runtimes** - Charts included with all batteries.
- + **Supports all PAP Devices** - No additional cables or kits.
- + **FAA Compliant** - Travel light and with ease.
- + **Compatible with all machines** - No additional purchases when upgrading machines.
- + **Designed with you in mind** - incorporates the latest AC power, USB & USB-c technology to power essential daily devices and gadgets.
- + **Versatile charging** - All batteries include a Zopec fast car charger and adapter to easily power up & go. Solar panels 40W & 60W are also available (not included).
- + **Higher capacity models available** - for any therapy configuration including; bi-level machines, humidifier and heated tubing use.
- + **Unlimited linking technology** - Connect as many batteries of the same model to increase capacity.



10044 Flanders CT NE, Suite 300, Blaine, MN 55449, USA.
Email: info@zopec.com | Phone: 612-234-7806
www.zopec.com



ZOPEC MEDICAL™
Innovation for Health

Sleep Safe. Sleep Uninterrupted.

Specifications	Explore 4000	Explore 5700	Explore 8000
Model #	214000C	215700C	218200C
Model Type	CPAP/BiPAP Only 1-2 nights	CPAP/BiPAP Only 2-3 nights	CPAP/BiPAP with humidifier (8-12 hrs.) or (24-32 hrs.) without
Power Outage Auto Detection	Yes	Yes	Yes
Battery Capacity		(160Wh)	(240Wh)
Dimensions Inches (L x W x H)	7.5 x 5.5 x 1	8.5 x 6 x 1	10 x 6 x 3
Weight (lbs.)	2	2.5	4
AC Output & Number of AC Outlets	120W	120W	250W
	1	1	2
USB-C	15 Watts (3A)	15 Watts (3A)	15 Watts (3A)
USB-A	15.5 Watts (3.1A)	15.5 Watts (3.1A)	15.5 Watts (3.1A)
Air Travel	Yes	Yes	No
Recharge Time	2-3 Hours	4-5 hours	7-8 hours
HSA/FSA Eligible	Yes	Yes	Yes
Charge Indicator	Digital Display	Digital Display	Digital Display
Battery Life	500/2500	500/2500	500/2500

All Explore batteries include:

- + Zopec fast car charger
- + Power wall adapter
- + 1 year warranty



Higher capacity models available. Ask us for details.
Available to ship globally.

Email: info@zopec.com | Phone: 612-234-7806
www.zopec.com