


### Fitting Guide for AirFit P10



**Note:** For complete user instructions, refer to the AirFit P10 User Guide. Rx Only

\*Declared dual numbered noise emission values (with uncertainty, 3dB(A) in accordance with ISO 8871.  
1. ResMed Testing Ref 0000-0201; 4.80298-206. 2. ResMed Clinical Study comparing AirFit P10 to Swift FX Ref 0629-117. ResMed objective data (pressure, leak, usage, AHI).  
ResMed Corp San Diego, CA USA +1 858 838 5000 or 1 800 424 0737 (toll free). ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 858 189 (toll free). See ResMed.com for other ResMed locations worldwide. AirFit, QuickFit and QuietAir are trademarks of ResMed Ltd. Swift is a trademark of ResMed Ltd and is registered in the U.S. Patent and Trademark Office. ©2013 ResMed. Specifications may change without notice. 1016801/3 2013-12

ResMed.com



## A quieter mask to fit your sleep.

The AirFit™ P10 nasal pillows mask is part of ResMed's new AirFit series—the right fit for you. This is our quietest mask yet, offering the reliability and performance of the Swift™ FX in a design that's approximately 50% lighter and 50% quieter.\* A study also found that patients using the AirFit P10 slept more than 40 minutes longer per night.<sup>2</sup>

#### Quieter

QuietAir™ woven-mesh vent technology makes this ResMed's quietest mask yet, promoting patient acceptance.

#### Lighter

50% lighter than our most popular nasal pillows mask,<sup>1</sup> offering patients the comfort of a better night's sleep.

#### Easier to Use

Key-clip nasal pillows feature color-coded sizing and left/right orientation for fast and easy fitting.

ResMed.com

## More than 40 minutes of extra sleep, every night

In a recent clinical study, we found that patients slept more than 40 minutes longer per night when using the ResMed AirFit P10.<sup>2</sup> That's almost five hours per week of additional sleep, helping ensure your patients will stay compliant with their sleep therapy.



### AirFit P10 Features and Benefits

**Dual-wall technology** — with trampoline action seals on contact, dispersing air gently to the nose.


**Key-clip nasal pillows** — offer effortless assembly, with color coding for easy sizing and left/right orientation to ensure patients fit the mask properly.

**Innovative flexible tube** — is ultra light and flexible yet highly durable, offering patients full freedom of movement without compromising seal.

**QuietAir woven-mesh venting** quietly directs exhaled air, so it's barely noticeable to patients and their bed partners.

**The AirFit P10 for Her** is designed especially for women, offering a smaller size designed for female facial contours, along with stylish tones and highlights.

### AirFit P10 Components



**3-part simplicity** makes the AirFit P10 50% lighter than the leading nasal pillows mask.<sup>1</sup>

**XS** Extra Small  
**S** Small  
**M** Medium  
**L** Large

#### Product Codes

AirFit P10 - complete system*	<b>62900</b>
AirFit P10 for Her - complete system*	<b>62910</b>

\*S, M and L pillows come packaged with the mask.  
\*\*XS, S and M pillows come packaged with the mask.

Note: For details, refer to the AirFit P10 User Guide.

\*Combined responses of those expressing a preference in an independent patient survey in the USA, UK, Germany, France & Australia, from September 2012. Visit [www.resmed.com/maskbrand](#).

