

ResMed

Actual Size: 140 mm (Width)
x 175 mm (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size


AirFit™ N30i

Nasal cradle mask



© 2018 ResMed Ltd. All Rights Reserved.
Specifications may change without notice.
Product availability may vary across regions.
10111103/1 2018-12

Rx Only CE 0123



Sample size check

10 mm 20 mm 30 mm 40 mm 50 mm 60 mm 70 mm 80 mm 90 mm 100 mm 110 mm 120 mm 130 mm 140 mm 150 mm 160 mm 170 mm