





**Ultra Mirage™ II**  
NASAL MASK


**Mirage Activa™**  
NASAL MASK




↑  
**STANDARD**



↑  
**LARGE**



↑  
**SHALLOW**



↑  
**SHALLOW WIDE**  
(Ultra Mirage II only)

**Actual Size: 5 1/2" (Width) x 7" (Height)**  
Ruler is for scaling purposes.



**CPAPMAN.COM**  
BEST SERVICE SINCE 1997

**Click Here To Order Today**

How to use the fitting template:

1. Cut out the template along the dotted line
2. Cut out nasal holes  
*(hint: fold each nasal hole in half to cut)*
3. Follow directions on template to determine appropriate cushion size

**RESMED**