

PHILIPS

RESPIRONICS

Sleep therapy

PAP retail accessories



CPAPMAN.COM
BEST SERVICE SINCE 1997

Click Here To Order Today

Rediscover dreams



A man in a blue sweater is shown from the chest up, looking down at a white PAP machine he is placing into a blue travel case. The case is open on a bed with a white blanket and a brown crocheted blanket. The background is a softly lit room with a lamp and a window.

PAP accessories

to meet your patients' needs

Discover new ways to enhance your patients' therapy experience by offering a wide variety of convenient PAP accessories. Carrying these accessories can enable you to:

- Keep in contact with patients and build lasting relationships for future opportunities
- Help patients enjoy a better quality of life with accessories that allow them to conveniently travel with their sleep therapy system and easily and discreetly store it at home
- Create a retail revenue opportunity and attract additional patients to your business

DreamStation travel case

The DreamStation travel case* gives patients a convenient, compact, and discreet way to travel with their DreamStation sleep therapy system and supplies.

Convenient




The case includes a DreamStation 65W power supply that is lighter and smaller than the DreamStation 80W power supply. The included USB port powers the DreamStation** system and can be used to charge items such as a smartphone and other electronic devices during travel.

Compact

The lightweight, hard-fabric case offers an ideal combination of protection and portability. Its individual compartments accommodate the DreamStation** system and neatly organizes the supplies. An included mesh pouch holds the tubing securely in place.

Discreet

The case is designed to make traveling with sleep therapy equipment less cumbersome. A molded insert liner helps the patient thoroughly pack the PAP system components. The case can be easily attached to a roller suitcase and its zip top helps facilitate airport screening. On a plane, it easily fits in the overhead bin or under the seat.



Benefits for patients

- The hard-fabric case provides added peace of mind that the sleep therapy system will be protected during travel.
- The included power supply offers travel convenience.
- The USB port on the 65W power supply reduces the need for an additional charger.
- Molded insert liner designed to help the patient thoroughly pack the PAP device components.

Components

The case accommodates the Philips Respironics DreamStation sleep therapy system (blower only, no humidifier), the patient's mask and tubing, a 65W power supply w/ USB port, and a removable luggage strap.

Specifications

Dimensions: 10.42" (26.46 cm) x 3.94" (10 cm) x 13.78" (35 cm)

Ordering information

Part number

Travel case, 65W power supply w/ USB port	1120135
Travel case, 65W power supply w/ USB port, 5/case	1124301

*DreamStation, mask, and tubing not included

**Without humidifier



PAP travel briefcase with detachable laptop bag

The PAP travel briefcase* gives patients an all-in-one carrying option that accommodates most PAP therapy systems and humidifiers, a laptop, and other small items. It makes it easy and efficient to travel with a sleep therapy system.

The travel briefcase features two bags that conveniently zip together into one carry-on. The therapy system and accessories fit snugly and securely into customizable compartments on one side of the case, while the laptop and its accessories fit securely into the ample-sized pockets of the other side. There's even space for books, magazines, and other small travel items.

Benefits for patients

- Its compartments protect the sleep therapy system and accessories during travel.
- Its zip top facilitates easy airport screening process.
- It fits within the FAA maximum carry-on size.
- The lightweight case offers easy portability.
- It's designed to easily attach to a roller suitcase.



Power supply

Tubing

DreamStation
with
humidifier



Removable
laptop bag

Detachable
shoulder strap

Components

The briefcase accommodates most PAP therapy devices and accessories and standard-size laptop computers.

Specifications

Dimensions: 16.14" w (40.99 cm) x 13.58" h (34.49 cm) x 7.87" d (19.99 cm)
Color: Black

Ordering information

PAP travel briefcase

Part number

1114784

*DreamStation, power supply and tubing not included



Mask and tubing bedside organizer

The bedside organizer gives patients a way to conveniently store their mask and tubing without cluttering a nightstand – both at home and during travel. It also keeps the components dust-free and ready for use.

The organizer features a weighted flap to allow it to hang over a nightstand or between a patient's mattress and box spring. It even features a zippered compartment for the mask and tubing, making the organizer a great travel companion.

Benefits for patients

- The organizer provides a clean, discreet place to store mask and tubing while not in use.
- Its design looks right at home in the bedroom.
- It can be used on a nightstand or in between mattress and box spring.

Specifications

Dimensions: 7.99" (20.3 cm) x 10.94" (27.8 cm) x 3.74" (9.5 cm)

Color: Light gray

Ordering information

Part number

Bedside organizer

1122135

Bedside organizer, 10/case

1124224



LiquiCell® nasal CPAP cushions

LiquiCell® is a thin, pliable cushion infused with fluid that can be applied virtually anywhere on the face, making it versatile. The cushion enhances comfort by acting as a barrier between the user's face and the mask. Gentle adhesive keeps the cushion securely in place.

Ordering information

Part number

LiquiCell nasal CPAP cushions, pack of 30

1109097

*LiquiCell is a registered trademark of Birchwood Laboratories, Inc.



Order your PAP accessories today

For more information, call 1-800-345-6443.



© 2016 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.com/respironics

Caution: US federal law restricts these devices to
sale by or on the order of a physician.

edoc EB 3/6/16 MCI 4107205
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 • 724-387-4000